


## PAHs

### Polyaromatic Hydrocarbons

CHROM  
9400




- |                   |                             |
|-------------------|-----------------------------|
| 1. Benzene        | 10. Benz[a]anthracene       |
| 2. Naphthalene    | 11. Chrysene                |
| 3. Acenaphthylene | 12. Benzo[b]fluoranthene    |
| 4. Fluorene       | 13. Benzo[k]fluoranthene    |
| 5. Acenaphthene   | 14. Benzo[a]pyrene          |
| 6. Phenanthrene   | 15. Dibenzo[a,h]anthracene  |
| 7. Anthracene     | 16. Benzo[g,h,i]perylene    |
| 8. Fluoranthene   | 17. Indeno[1,2,3,c,d]pyrene |
| 9. Pyrene         |                             |

**Column:** Alltech® Prevail™ C18, 5µm, 150 x 4.6mm (Part No. 99208)  
**Mobile Phase:** A: Water B: Acetonitrile  
**Gradient:**

Time:	0	4	15
%B:	72	72	100

  
**Flow Rate:** 1.5mL/min  
**Detector:** UV at 254

### Beyond the EPA 16 Priority-Pollutant PAHs



- |                        |                                    |                             |
|------------------------|------------------------------------|-----------------------------|
| 1. Naphthalene         | 9. Fluoranthene                    | 17. Benzo[e]pyrene          |
| 2. Acenaphthylene      | 10. Pyrene                         | 18. Benzo[b]fluoranthene    |
| 3. 1-Methylnaphthalene | 11. Benzo[c]phenanthrene           | 19. Benzo[k]fluoranthene    |
| 4. 2-Methylnaphthalene | 12. Cyclopenta[c,d]pyrene          | 20. Benz[a]pyrene           |
| 5. Acenaphthene        | 13. Benz[a]anthracene              | 21. Dibenzo[a,h]anthracene  |
| 6. Fluorene            | 14. Chrysene                       | 22. Benzo[g,h,i]perylene    |
| 7. Phenanthrene        | 15. Benzo[b]naphtho[2,1-d]thiopen  | 23. Indeno[1,2,3,c,d]pyrene |
| 8. Anthracene          | 16. 7,12-Dimethylbenz[a]anthracene |                             |

**Column:** Vydac® C18, 5µm, 150 x 4.6mm (Part No. 201TP5415)  
**Mobile Phase:** A: Water B: Acetonitrile  
**Gradient:**

Time:	0	30
%B:	50	100

  
**Flow Rate:** 1.0mL/min

### Fast Polyaromatic Hydrocarbons

CHROM  
10757

- |                   |                          |                             |
|-------------------|--------------------------|-----------------------------|
| 1. Naphthalene    | 7. Fluoranthene          | 13. Benzo[a]pyrene          |
| 2. Acenaphthylene | 8. Pyrene                | 14. Dibenzo[a,h]anthracene  |
| 3. Acenaphthene   | 9. Benzo[a]anthracene    | 15. Benzo[g,h,i]perylene    |
| 4. Fluorene       | 10. Chrysene             | 16. Indeno[1,2,3,c,d]pyrene |
| 5. Phenanthrene   | 11. Benzo[b]fluoranthene |                             |
| 6. Anthracene     | 12. Benzo[k]fluoranthene |                             |



**Column:** Alltech® Alltima™ HP C18 HiLoad, 3µm, 53 x 7mm (Part No. 87690)  
**Mobile Phase:** A: Water B: Acetonitrile  
**Gradient:**


Time:	0	0.1	5
%B:	40	40	100

  
**Flow Rate:** 5mL/min  
**Detector:** UV at 230nm

### Nitroaromatics (Normal-Phase)

CHROM  
10358


- |                       |
|-----------------------|
| 1. Nitrobenzene       |
| 2. 2,6-Dinitrotoluene |
| 3. 1,3-Dinitrobenzene |
| 4. 4-Nitrophenol      |



**Column:** Alltech® Alltima™ HP CN, 5µm, 150 x 4.6mm (Part No. 87769)  
**Mobile Phase:** Hexane:Isopropanol (85:15)  
**Flow Rate:** 1.0mL/min  
**Detector:** UV at 254nm

# Polycyclic aromatic hydrocarbons (PAH's)


## 10 036 Isocratic Separation of Polycyclic Aromatic Hydrocarbons (PAH's)



**Column phase:** GROM-SIL PAH -1 (tailor-made)  
**Column size:** 250 x 4mm  
**Eluent:** water / ACN = 15 / 85  
**Flow rate:** 1.5 ml/min  
**Pressure:** 13 MPa  
**Temperature:** 25°C  
**Detection (UV):** 254 nm  
**Injection:** 5 µl (5 µg/ml each)

- |                   |                            |
|-------------------|----------------------------|
| 1) Naphthalene    | 9) Benz(a)anthracene       |
| 2) Acenaphthylene | 10) Chrysene               |
| 3) Acenaphthene   | 11) Benzo(b)fluoranthene   |
| 4) Fluorene       | 12) Benzo(k)fluoranthene   |
| 5) Phenanthrene   | 13) Benzo(a)pyrene         |
| 6) Anthracene     | 14) Dibenzo(a,h)anthracene |
| 7) Fluoranthene   | 15) Benzo(g,h)perylene     |
| 8) Pyrene         | 16) Ideno(1,2,3,c,d)pyrene |


## 10 163 Analysis of Polycyclic Aromatic Hydrocarbons (EPA Standard)



**Column phase:** GROM-PAH, 3 µm  
**Column size:** 50 x 4 mm  
**Eluent A:** H<sub>2</sub>O / MeOH = 20 / 80  
**B:** ACN  
**Gradient:** 10% B (0-0.5 min),  
 100% B (0-1.0 min),  
 100% B (1.0-4.0 min)  
**Flow rate:** 2.0 ml/min  
**Pressure:** 16.4 MPa  
**Temperature:** RT  
**Detection (UV):** 254 nm (1.3 µl flow cell)  
**Injection:** 2 µl (10 mg/ml of each)

- |                            |
|----------------------------|
| 1) Naphthalene             |
| 2) Acenaphthylene          |
| 3) Acenaphthene            |
| 4) Fluorene                |
| 5) Phenanthrene            |
| 6) Anthracene              |
| 7) Fluoranthene            |
| 8) Pyrene                  |
| 9) Benz(a)anthracene       |
| 10) Chrysene               |
| 11) Benzo(b)fluoranthene   |
| 12) Benzo(k)fluoranthene   |
| 13) Benzo(a)pyrene         |
| 14) Dibenzo(a,h)anthracene |
| 15) Benzo(g,h)perylene     |
| 16) Ideno(1,2,3,c,d)pyrene |

## 10 037 Separation of PAH's by Gradient Elution



**Column phase:** GROM-SIL PAH -1 (tailor-made)  
**Column size:** 250 x 4mm  
**Eluent A:** water  
**B:** ACN  
**Gradient:** 50% B (0-5 min), 50-85% B (5-15 min),  
 85-100% B (15-20 min),  
 100% B (20-25 min)  
**Flow rate:** 1.5 ml/min  
**Pressure:** 24 MPa  
**Temperature:** 25°C  
**Detection (UV):** 254 nm  
**Injection:** 5 µl (5 µg/ml, each)

- |                   |                            |
|-------------------|----------------------------|
| 1) Naphthalene    | 9) Benzo(a)anthracene      |
| 2) Acenaphthylene | 10) Chrysene               |
| 3) Acenaphthene   | 11) Benzo(b)fluoranthene   |
| 4) Fluorene       | 12) Benzo(k)fluoranthene   |
| 5) Phenanthrene   | 13) Benzo(a)pyrene         |
| 6) Anthracene     | 14) Dibenzo(a,h)anthracene |
| 7) Fluoranthene   | 15) Benzo(g,h)perylene     |
| 8) Pyrene         | 16) Ideno(1,2,3,c,d)pyrene |

Columns	Dimensions	Field of application	Order number
PAH-column for EPA-standard - (separation of the 16 PAH's)	250 x 2 mm cartridge 250 x 4 mm cartridge 50 x 4.6 mm cartridge	environmental analysis: soil, pollution (air / water), etc.	GS PA 1 0530 K 2502 GS PA 1 0530 K 2504 GS PA 1 0330 K 0505

**Note!** All products are also available as columns (complete with end-fittings). Please enquire for other dimensions. Reduced prices for Refill-columns.