

## VisionHT™ Columns

### One Platform for UHPLC, HPLC and PREP

- Easy method transfer between UHPLC, HPLC and preparative chromatography
- Broad selectivity range—unique complementary phase chemistries
- Maximize speed and resolution with sub 2µm on any LC system





Separate complex samples with excellent resolution with VisionHT™ high performance columns. They are available with 1.5, 3, 5 and 10µm particle sizes and six different phase selectivities spanning the full polarity spectrum. The same base silica and bonding chemistries are used across all particle sizes therefore offering simple method transfer between UHPLC, HPLC and preparative systems.

VisionHT™ Phase Specifications							
Phase	Pore Size (Å)	Surface Area (m <sup>2</sup> /g)	Carbon Load (%)	Endcapped	Chromatographic Properties	Application/Benefit	USP L-code
C18 HighLoad	120	220	11	Yes	Ultra-high purity silica, fully bonded.	General purpose for broad range compounds, classic selectivity, high-capacity for hydrophobic compounds.	L1
C18 Basic	120	220	5	Proprietary	Ultra-high purity silica. Controlled silica surface exposure gives dual mode separation with polar and non-polar analytes.	Alternate reversed-phase selectivity. High polar retention especially with compounds having two or more polar groups. Excellent sensitivity and peak shape for basic compounds, without the need for acidified mobile phases.	L1
C18 Classic	100	200	6	Yes	Lower carbon load. Slight silica exposure.	Reversed-phase separations with reduced bonding optimized for speed. Some additional polar retention.	L1
C18 Polar	100	200	5	No	High silica exposure, low carbon load. Uniform coverage of inert vicinal silanols.	Unique polar selectivity. Low carbon load gives fastest reversed-phase elution times while retaining polar compounds longer.	L1
HILIC	120	220	—	No	Polar phase with shorter equilibration times. Shipped in ACN/Water.	Peak reversal compared to reversed-phase. Ideal for very polar compounds with high organic mobile phases for improved sensitivity by MS.	L3
Silica	120	220	—	No	Traditional normal phase for use in 100% organic mobile phases.	For isomeric separation of non-aqueous compatible compounds by absorption chromatography.	L3

### Simple Method Transfer




Optimizing and transferring methods between HPLC, UHPLC and preparative systems is not a simple and intuitive task. However, when the same base silica and bonding is available in sub-2, 3, 5, and 10µm particles and in a variety of column formats, methods can be more easily transferred across system types and between laboratories, improving efficiency and productivity.

### Identical Selectivity Across All Particle Sizes



## Broad Selectivity Range

Four reversed-phase column chemistries offer complementary selectivity. The unique VisionHT™ C18 Polar and Basic phases retain polar compounds while the C18 HighLoad has a more traditional reversed-phase selectivity. Using phases with orthogonal selectivity gives confidence that if one phase does not produce the separation the other will.




# hplc columns & accessories


hplc columns

## Consistent and Reliable Performance


Ultra-high mechanical strength silica and robust bonding chemistry minimize variations in capacity and selectivity. Media consistency combined with reproducible column packing methods delivers reliable performance and long column lifetime.




VisionHT™ silica's surface area, pore volume, and pore diameter are highly reproducible.



VisionHT™ media and advance packing methods produce columns with highly consistent capacity and selectivity.



Before and after chromatograms show a consistent level of performance after exposure to 16,000psig pressures for 20,000 column volumes.



Competitive columns lose performance over time under high-pressure conditions while VisionHT™ columns remain stable.

## Method Development Ultra MD® Kits — Experience Broad Selectivity Range of the VisionHT™ Phases

Ultra MD® kits include 4 phases with orthogonal selectivity to enable fast and efficient method development for the separation of very polar to very hydrophobic compounds.

### VisionHT™ Column Ultra MD® Kits

Description	Phases	Dimensions	Part No.
Ultra MD® Kit 1	C18 Classic, C18 Polar, C18 HighLoad, C18 Basic	2 x 100mm	5142692
Ultra MD® Kit 2	C18 Classic, C18 Polar, C18 HighLoad, C18 Basic	2 x 50mm	5142693
Ultra MD® Kit 3	C18 Classic, C18 Polar, C18 HighLoad, C18 Basic	1 x 50mm	5142691

## VisionHT™ Guard Columns — Protect Your Column to Minimize Downtime and Reduce Cost



Couple the zero-dead volume integral guard directly to any VisionHT™ column, with no loss in efficiency.

hplc columns

### VisionHT™ Columns

Format		i.d. (mm)	Length (mm)	Particle Size (µm)	C18 HighLoad Part No.	C18 Basic Part No.	C18 Classic Part No.	C18 Polar Part No.	HILIC Part No.	Silica Part No.
UHPLC	U	1	20	1.5	5142540	5141902	5139555	5139556	5141910	—
		1	30	1.5	5142541	5141903	5139559	5139601	5141912	—
		1	50	1.5	5142542	5141904	5139603	5139604	5141913	5141921
		1	100	1.5	5142543	5141905	5139607	5139608	5141914	5141923
		2	20	1.5	5142544	5141906	5139557	5139558	5141916	—
		2	30	1.5	5142545	5141907	5139600	5139602	5141917	—
		2	50	1.5	5142546	5141908	5139605	5139606	5141919	5141922
		2	100	1.5	5142547	5141909	5139609	5139610	5141920	5141924
		2	150	1.5	5148861	5148862	5148863	5148864	5149010	—
		Expedite™	E	2.1	20	1.5	5151986	5152009	5152043	5152077
4.6	20			1.5	5151985	5152008	5152042	5152076	—	—
Microbore	M	2.1	50	3	5151981	5152004	5151948	5152072	5152091	—
		2.1	100	3	5151980	5152003	5151947	5152071	—	—
		2.1	150	3	5151979	5152002	5151946	5152070	5152090	5152116
Solvent Reducer	SR	3	100	3	5151978	5152001	5152025	5152069	—	—
		3	150	3	5151977	5152000	5152024	5152068	—	—
		3	250	3	5151975	5151998	5152022	5152066	5152089	5152115
		3	150	5	5151976	5151999	5152023	5152067	—	—
Analytical	A	3	250	5	5151974	5151997	5152021	5152065	5152088	5152114
		4.6	50	3	5151973	5151996	5152020	5152064	5152087	5152113
		4.6	100	3	5151972	5151995	5152019	5152063	—	—
		4.6	150	3	5151971	5151994	5152018	5152062	5152086	5152112
		4.6	150	5	5151970	5151993	5152017	5152061	5152085	5152111
Rocket™	R	4.6	250	5	5151920	5151921	5152016	5152060	5152084	—
		7	33	1.5	5151984	5152007	5152041	5152075	5152110	—
		7	53	3	5151983	5152006	5152040	5152074	5152093	—
		7	53	1.5	5151982	5152005	5151949	5152073	5152092	—
Preparative	P	10	150	5	5151987	5152010	5152044	5152078	—	—
		10	250	5	5151988	5152011	5152045	5152079	—	—
		22	150	5	5151989	5152012	5152046	5152080	—	—
		22	250	5	5151990	5152013	5152047	5152081	—	—
		22	250	10	5151991	5152014	5152048	5152082	—	—
Capillary Guard*	G	1	5	1.5	5142549	5141953	5141950	5141951	5141955	5141957
All-Guard™**	G	2	5	1.5	5142548	5141952	5141594	5141595	5141954	5141956
		2.1	7.5	3	5155906	5155908	5155907	5155909	5155911	5155912
		4.6	7.5	5	5155900	5155902	5155901	5155903	5155904	5155905

\*Requires Capillary holder, Part No. 3118351.

\*\*Requires All-Guard™ Holder, Part No. 80101.

### more info

For detailed information on column formats, see page 3.

### technical assistance

Contact Tech Support: Email: [discoverysciences@grace.com](mailto:discoverysciences@grace.com)

Online: [www.discoverysciences.com](http://www.discoverysciences.com)