

## Acceptable Modifications to an HPLC Validated Method

Even if you are using an FDA validated or a USP recommended method, some operating conditions can be adjusted if the modifications respect the acceptable specifications proposed by Pharmacopeias<sup>1-3</sup> and the FDA<sup>4</sup>. A side-by-side comparison of both the original and the adjusted method needs to be performed to demonstrate that the method's accuracy and precision is not affected by these modifications.

Acceptable Modifications to an HPLC Validated Method		
Parameters	Allowable modification	Examples of possible modifications
Mobile phase pH	± 0.2 units	Validated pH: 7.0 Allowed pH range: 6.8 - 7.2
Concentration of salts in buffer	± 10%	Validated concentration: 20 mM Allowed concentration range: 18 - 22 mM
Ratio of components in mobile phase	Only the minor components can be adjusted by ± 30% or ± 2% absolute ( <i>i.e.</i> : in regards to the total mobile phase), whichever is the larger but should never exceed ± 10% absolute or removed totally.	<b>Binary mixtures:</b> Validated ratio: 50/50 Allowed ratio: 40/60 to 60/40  Validated ratio: 95/5 Allowed ratio: 93.5/6.5 to 96.5/3.5  <b>Ternary mixtures:</b> Validated ratio: 60/35/5 Allowed % of the 1 <sup>st</sup> component: 60% Allowed % of the 2 <sup>nd</sup> component: 25 - 45% Allowed % of the 3 <sup>rd</sup> component: 3.5 - 6.5%  The total of the three components together need to be 100%.
Wavelength of UV detector	No modification allowed.	n/a
Column length	± 70%	Validated length: 150 mm Allowed length range: 45 - 255 mm
Column inner diameter	± 50%	Validated inner diameter: 4.6 mm Allowed inner diameter range: 2.3 - 10.6 mm
Flow rate	± 50%	Validated flow rate: 1.00 mL/min Allowed flow rate range: 0.5 - 1.5 mL/min
Injection volume	May be increased to as much as 2 times if no adverse effects on LOD and repeatability.	n/a
Particle size	No increase permitted. May be decreased by as much as 50%.	Validated particle size: 5 µm Allowed particle size range: 2.5 - 5 µm
Column temperature	± 20%	Validated temperature: 23°C Allowed length range: 18.4 - 27.6°C

<sup>1</sup> USP. USP 32-NF 27, Chromatography <621>. Rockville, MD: USP; 2009:227.

<sup>2</sup> USP. Second Supplement to USP 32-NF 27. Rockville, MD: USP; 2009:4147.

<sup>3</sup> USP. USP 32-NF 27, Verification of Compendial Procedures <1226>. Rockville, MD: USP; 2009:736.

<sup>4</sup> ORA Laboratory Procedure, Food and Drug Administration, modification criteria.