

AN number: LA1011
Date: 2015-04-16
Revision: 1.0

The determination of nitrogen in pet food

As the Kjeldahl method involves use of highly corrosive and hazardous chemicals at high temperature it is strongly recommended to use protective glasses and gloves. The safety sections in the KjelROC Operation Manuals should be read before starting any analysis work.

INTRODUCTION

Protein determination is one of the key analyses performed in the food and feed industry. The samples require digestion with sulfuric acid to convert nitrogen to ammonium sulfate. An excess of sodium hydroxide is added to the digest to liberate ammonia. The sample is steam distilled and the liberated ammonia is collected in a boric acid solution and titrated with hydrochloric acid, using a colorimetric end-point detection. The measured nitrogen content is calculated from the amount of determined ammonia.

EXPERIMENTAL

Apparatus

- KjelROC Digester: DI-210-A, DI-220-A, DI-310-A or DI-320-A
- Exhaust system, compatible with the Digester
- KjelROC Analyzer KD-310-A (steam distillation/auto titration) or KjelROC Distillation Unit KD-210-A (steam distillation/manual titration)
- KjelROC Scrubber DI-110-A or KjelROC Digester-water jet pump DI-054-A
- Analytical balance
- Digestion tubes – 250 ml

Reagents and accessories

- Sulfuric acid, concentrated (95-98%)
- Kjeldahl tablets, Missouri Tablet, Cu 0,3%, 3,5 g: KT-211-A
- Sodium Hydroxide, 40% (w/w)
- Hydrochloric acid, 0,2000 N
- Boric acid, 1% (w/v). See LA1000 for preparation

Sample preparation

Prepare the sample by using a suitable mixer. Weigh 1 g of the homogenized sample into a test tube.

Digestion

Add 2 Kjeldahl tablets and 10 ml H₂SO₄ to digestion tube. Shake gently to “wet” the sample. Place tubes in digestion block (420°C) and position the exhaust and turn on the Scrubber or water jet pump. Digest for 60 minutes. Remove rack with exhaust and let cool.

Distillation

On some systems part or all of this is performed automatically.

Dilute cooled digest with 70 ml H₂O and add 50 ml NaOH (40%). Add 30 ml boric acid to receiver vessel. Distil and titrate with HCl (0,2000 N). Perform a blank before each batch of samples.

CALCULATIONS

$$\% \text{ Nitrogen} = \frac{(T - B) \times N \times 14,007 \times 100}{\text{weight}_{\text{sample}}(\text{mg})}$$

T = Sample titration

B = Blank titration

N = Normality of titrant

$$\% \text{ Protein} = N \times F$$

F = protein factor = 6,25 (5,70 for wheat)

REFERENCES

This Application Note should be used in conjunction with Application Note LA1000 “Application Guide”

OG 1000 KjelROC Analyzer Operation Manual

OG 1010 KjelROC Distillation Unit Operation Manual

OG 1004 KjelROC Digestor Operation Manual

OG 1006 KjelROC Scrubber Operation Manual

