



Solutions for Feed

Wet Chemistry Analyses for Producers of Animal and Fish Feed

Crude Protein

OP SIS LiquidLINE has solutions for determination of Kjeldahl (TKN) protein following standard methods.

The samples are digested with sulphuric acid to convert nitrogen into ammonium sulphate. The samples are further distilled by steam distillation followed by titration.

Examples: Protein in animal feed stuff, forage (plant tissue), grain, oilseeds, pet food and protein in residues from other production processes such as palm kernel cake.

Our Solution

- The KjelROC Digestor Advanced motor lift makes the digestion efficient and saves valuable operator time.
- KjelROC Analyzer with integrated Titration offers titration with low relative standard deviation and wireless communication saving time and costs.

Standards

ISO 1871:2009, ISO 5983-2:2009
AOAC 2011.11
AOAC 954.01
AOAC 976.06
AOAC 984.13

Application Notes

LA1000 Application Guide Kjeldahl
LA1011 Determ. of nitrogen in pet food
Further Notes on request

Total Fat

OP SIS LiquidLINE provides instruments to determine Total Fat according to standard methods.

The sample is hydrolysed and thereafter extracted in hot solvents. Calculation of total fat content follows after the extract has been dried to a constant weight.

Examples: Fat in Cattle, Cat, Dog, Horse, Chicken, Rabbit and Pig feed

Our Solution

- The HydROC hydrolysis unit offers a unique filter technology that saves time and reduces the risk of errors when moving samples between hydrolysis and extraction.
- The SoxROC extraction unit with batch handling and full automation facilitates the extraction.

Standards

98/64/EC
ISO 6492
ISO/FDIS 11085

Application Notes

LA1002, Appl. Guide Solvent Extraction
LA1005, Extr. of Total Fat in Petfood
Further Notes on request

Crude Fat

OP SIS LiquidLINE provides instruments to determine Crude Fat with Hot Solvent extraction.

The sample is prepared and thereafter extracted in hot solvents. Calculation of fat content follows after the extract has been dried to a constant weight.

Examples: Fat in pet food, feeds, cereal grains and forages. Fat in fish meals.

Our Solution

- The SoxROC extraction unit with batch handling and full automation facilitates the extraction.
- The instrument provides significant time savings versus cold extraction and a recovery of over 90% of used solvents.

Standards

AOAC 920.39
AOAC 948.15, AOAC 948.16
AOAC 954.02
AOAC 2003.05, AOAC 2003.06
ISO/FDIS 11085

Application Notes

LA1002, Appl. Guide Solvent Extraction
Further Notes on request

Oil from distilled grains waste

When producing ethanol there are waste products from the distilled grains. The waste from this process is typically sold as cattle feed whereas the price is based on the oil content of the waste. OPSIS LiquidLINE provides instruments to extract oil from this waste.

Our Solution



- The SoxROC extraction unit with batch handling and full automation facilitates the extraction.
- The instrument provides significant time savings versus cold extraction and a recovery of over 90% of used solvents.

Standards

-

Application Notes

LA1002, Appl. Guide Solvent Extraction
Further Notes on request



Azerbaijan's leading livestock feed producer is using a KjelROC System in their laboratory. The firm's priority is to bring high productivity, modern standards and a developed structure to the largely traditional livestock farming in Azerbaijan.

Premix Factory No.1, Biotechnological Center, outside Belgorod, Russia is a factory that produces Lysine. Lysine production for animal feed is a major global industry and an important additive to animal feed. The laboratory is using a complete OPSIS LiquidLINE system with one KjelROC Analyzer, one KjelROC Digester 10 Advanced and a KjelROC Scrubber.



State Veterinary Institute Prague, Czech Republic, is specialized in veterinary and laboratory diagnostic of animal feed. They use a full KjelROC Analyzer system with Digester block and Scrubber.



OPSIS LIQUIDLINE - INNOVATIVE WET CHEMISTRY

OPSIS AB, founded in 1985 in Sweden, took the concept of measuring gases with light and developed it into a commercially viable product. In 2013, we took another step and moved our innovative technology into Wet Chemistry and Liquids. We can offer:

- AN APPLICATION LABORATORY READY TO ASSIST
- CUSTOMISED TRAINING AND SUPPORT FROM SWEDEN
- THE LATEST IN MAINTENANCE
- A COMPLETE PORTFOLIO

OPSIS AB, Box 244
SE-244 02 Furulund Sweden
Telephone +46 46 72 25 00
Telefax +46 46 72 25 01
E-mail info@opsis.se
www.liquidline.se

LB 1013
2015 05