


ANALYSIS OF AMINES ON BPX5

Column Part No.: 054113

Phase: BPX5, 0.25 μ m
Column: 25 m x 0.22 mm I.D.
Initial Temp.: 40 $^{\circ}$ C, 1 min
Rate: 8 $^{\circ}$ C/min
Final Temp.: 360 $^{\circ}$ C
Detector: FID, 380 $^{\circ}$ C
Injector: Split
Carrier Gas: H₂, 12 psi

Notes: BPX5 is a chemically inert column ideal for aromatic amine analysis




For more information contact our technical customer support team on: techsupport@sge.com

ANALYSIS OF AROMATIC AMINES ON BP5

Column Part No.: 054197

Phase: BP5, 1.0 µm film
Column: 12 m x 0.53 mm ID
Initial Temp: 60 °C, 0 min
Rate: 10 °C/min
Final Temp: 190 °C, 0 min
Detector: FID
Sensitivity : 128 x 10⁻¹² AFS
Injection Mode: Split




For more information contact our technical customer support team on: techsupport@sge.com

ANALYSIS OF AMINES ON BP1

Column Part No: 054097

Phase: BP1, 3.0 μm film
Column: 12 m x 0.53 mm ID
Initial Temp: 70 $^{\circ}\text{C}$
Rate: 10 $^{\circ}\text{C}/\text{min}$
Final Temp.: 250 $^{\circ}\text{C}$
Carrier Gas: Nitrogen
Injection Volume: 0.1 μL




For more information contact our technical customer support team on: techsupport@sge.com

ANALYSIS OF TRIETHYLAMINE AND TRIETHANOLAMINE ON SOLGEL-1MS™

Column Part No: 054798

Phase: SolGel-1ms™, 0.25 µm film
Amine mix: 10 ng/µl in dichloromethane
Column: 30 m x 0.32 mm ID
Initial Temp: 40 °C, 5.0 min
Rate 1: 20 °C/min to 200 °C
Final Temp: 200 °C, 7 min
Detector Type: FID
Detector Temp: 300 °C
Carrier Gas: Helium, 9.9 psi

Carrier Gas Flow: 2.2 mL/min
Constant Flow: On
Average
Linear Velocity: 35 cm/sec at 100 °C
Mode of Injection: Split
Split Ratio: 50:1
Injection Volume: 0.3 µL
Injection Temp: 250 °C
Autosampler: No



For more information contact our technical customer support team on: techsupport@sge.com

ANALYSIS OF ANILINE IMPURITIES ON BPX5


ANILINE IMPURITIES

Column Part No.: 054127

Phase: BPX5, 1.0 μm
Column: 12 m x 0.32 mm I.D.
Initial Temp.: 50 $^{\circ}\text{C}$, 1 min
Rate 1: 15 $^{\circ}\text{C}/\text{min}$
Final Temp.1: 120 $^{\circ}\text{C}$
Rate 2: 20 $^{\circ}\text{C}/\text{min}$
Final Temp.2: 330 $^{\circ}\text{C}$, 10 min

Detector: FID, 380 $^{\circ}\text{C}$
Injector: Split
Carrier Gas: He, 4 psi


Notes: Thick film BPX5 column provides excellent separation with high temperature compatibility



For more information contact our technical customer support team on: techsupport@sge.com

ANALYSIS OF AROMATIC AMINES ON BPX35

Column Part No:	Custom Made , for further information contact our technical support team	Constant Flow:	On
Phase:	BPX35 1.8 μm film	Average	
Amine mixture:	100 ng/ μL in dichloromethane	Linear Velocity:	41 cm/sec at 40 $^{\circ}\text{C}$
Column:	30 m x 0.32 mm ID	Injection Mode:	Split
Initial Temp:	40 $^{\circ}\text{C}$, 1 min	Purge on Time:	0.5 min
Rate 1:	20 $^{\circ}\text{C}/\text{min}$ to 240 $^{\circ}\text{C}$	Purge on (Split):	
Final Temp:	240 $^{\circ}\text{C}$, 7 min	Vent Flow:	60 mL/min
Detector Type:	Mass Spectrometer	Injection Volume:	0.1 μL
Carrier Gas:	He, 10.8 psi	Injection Temp:	250 $^{\circ}\text{C}$
Carrier Gas Flow:	1.3 mL/min	Liner Type:	4 mm ID Double Taper Liner
		Liner Part No:	092018




For more information contact our technical customer support team on: techsupport@sge.com

ANALYSIS OF AROMATIC AMINES FROM DIAZO DYES ON BPX35

Column Part No: 054701

Phase: BPX35 0.25 µm film
 Azo Dyes standard: 10 ppm solution in DCM
 Column: 30 m x 0.25 mm ID
 Initial Temp: 50 °C, 2 min
 Rate 1: 15 °C to 240 °C
 Rate 2: 10 °C to 280 °C
 Final Temp: 280 °C, 25 min
 Detector Type: MSD
 Carrier Gas: He, 7.1 psi
 Carrier Gas Flow: 1.0 mL/min
 Constant Flow: On

Average
 Linear Velocity: 36 cm/sec at 50 °C
 Injection Mode: Splitless
 Purge on Time: 1.0 min
 Purge on (Split)
 Vent Flow: 60 mL/min
 Injection Volume: 1 µL
 Injection Temp: 250 °C
 Liner Type: 4 mm ID Double Taper Liner
 Liner Part No: 092018



For more information contact our technical customer support team on: techsupport@sge.com