


# ANALYSIS OF $\alpha$ -IRONE (STANDARD) ON CYDEX-B

## $\alpha$ -IRONE (STANDARD)

<b>Column Part No.:</b>	<b>054900</b>	Rate:	2 °C/min
Column:	Cydex-B, 0.25 $\mu$ m	Final Temp.:	220 °C
Phase:	25 m x 0.22 mm I.D.	Carrier Gas:	Helium
Initial Temp.:	100 °C, 1 min	Carrier Pressure:	16 psi
		Injection Mode:	Split 50:1




For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

## ANALYSIS OF ORRIS CONCENTRATE ON BP20

## ORRIS CONCENTRATE

<b>Column Part No.:</b>	<b>054436</b>	Rate:	10 °C/min
Column:	BP20, 0.5 µm	Final Temp.:	250 °C, 10 min
Phase:	25 m x 0.32 mm I.D.	Carrier Gas:	Helium
Initial Temp.:	70 °C, 1 min	Carrier Pressure:	10 psi
		Injection Mode:	Split 50:1




For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# ANALYSIS OF LAVENDER OIL ON CYBEX-B

## LAVENDER OIL

<b>Column Part No.:</b>	<b>054901</b>	Detector:	FID
Phase:	Cydex-B 0.25 µm film	Sensitivity:	32 x 10 <sup>-12</sup> AFS
Column:	50 m x 0.22 mm ID	Injection Mode:	Split
Initial Temp.:	Isothermal at 90 °C		

Notes: The Cydex - B column enables easy detection of adulteration of Lavender oils




For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# ANALYSIS OF MENTHOL OIL ON CYDEX-B

## MENTHOL OIL

<b>Column Part No.:</b>	<b>054901</b>	Final Temp.:	130 °C
Phase:	Cydex-B, 0.25 µm film	Carrier Gas:	H <sub>2</sub>
Column:	50 m x 0.22 mm I.D.	Detector:	F.I.D.
Initial Temp.:	100 °C, 5 min	Sensitivity:	32 x 10 <sup>-12</sup> AFS
Rate:	2 °C/min	Injection Mode:	Split

Notes: *Cydex - B* column enables the separation of three different enantiomer pairs in Menthol Oil.




For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

## ANALYSIS OF BORNEOLS ON CYDEX-B

## ANALYSIS OF BORNEOLS


<b>Column Part No.:</b>	<b>054900</b>	Final Temp.:	130 °C
Phase:	Cydex-B, 0.25 µm film	Carrier Gas:	H <sub>2</sub> , 12 psi
Column:	25 m x 0.22 mm I.D.	Detector:	F.I.D.
Initial Temp.:	90 °C, 1 min	Sensitivity:	32 x 10 <sup>-12</sup> AFS
Rate:	2 °C/min	Injection Mode:	Split



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

## ANALYSIS OF RASPBERRY FLAVOURS ON BPX5


<b>Column Part No.:</b>	<b>054114</b>	Final Temp:	250 °C, 10 min
Phase:	BPX5, 0.25 µm film	Injector Temp.:	240 °C
Column:	50 m x 0.22 mm ID	Injector Mode:	Split
Initial Temp.:	60 °C, 10 min	Detector:	FID 320 °C
Rate:	2 °C/min	Carrier Gas:	Helium, 25 psi



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

## RASPBERRY EXTRACT ON BP1

<b>Column Part No.:</b>	<b>054066</b>	Temp 2:	130 °C
Phase:	BP1, 0.5 µm film	Rate 2:	4 °C/min
Column:	50 m x 0.32 mm ID	Final Temp.:	130 °C
Initial Temp.:	40 °C	Detector:	FID
Rate 1:	2 °C/min	Injection Mode:	Splitless (1 µL)




For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# ANALYSIS OF ROSEMARY ESSENTIAL OIL ON BPX5

**Column Part No.:** 054101  
**Phase:** BPX5, 0.25 µm film  
**Column:** 30 m x 0.25 mm ID  
**Initial Temp.:** 45 °C  
**Initial Hold Time:** 1 min.  
**Rate 1:** 5 °C/min  
**Temp 2:** 250 °C  
**Temp 1 Hold Time:** 15 min

**Detector:** 5973 MSD  
**Injector Mode:** Split (50:1)  
**Injector Temp.:** 250 °C  
**Carrier Gas:** He, 35 cm/sec at 45 °C  
**Injection Volume:** 1 µL  
**Concentration (ng/mL):** 1:20 dilution in DCM  
**GC model:** HP 6890 GC, 5973 MSD  
**Interface/Source/Quad:** 320 °C, 250 °C, 106 °C




For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)



# ANALYSIS OF EUCALYPTUS ESSENTIAL OIL ON BPX5

<b>Column Part No.:</b>	<b>054101</b>	Detector:	5973 MSD
Phase:	BPX5, 0.25 µm film	Injector Mode:	Split (50:1)
Column:	30 m x 0.25 mm ID	Injector Temp.:	250 °C
Initial Temp.:	45 °C	Carrier Gas:	He, 35 cm/sec at 45 °C
Initial Hold Time :	1 min	Injection Volume :	1 mL
Rate 1:	5 °C/min	Concentration (ng/mL)	1:20 dilution in DCM
Temp 2:	250 °C	GC model:	HP 6890 GC, 5973 MSD
Temp 1 Hold Time:	15 min	Interface/Source/Quad:	320 °C, 250 °C, 106 °C


*Note: The non-polar nature of the phase will generally give separation based on boiling point and this is complimentary to the unique selectivity of the polyethylene glycol BP20 phase*



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

## ANALYSIS OF NUTMEG OIL ON BPX5

<b>Column Part No.:</b>	<b>054101</b>	Average Linear Velocity:	36 cm/sec at 40 °C
Phase:	BPX5, 0.25 µm film	Injection Mode:	Split
Column:	30 m x 0.25 mm ID	Split Ratio:	200:1
Initial Temp.:	40 °C, 1 min.	Purge on (Split) Vent Flow:	200 mL/min.
Rate:	5 °C/min to 260 °C,	Injection Volume:	0.2 µL
Final Temp:	260 °C	Injection Temp.:	250 °C
Detector Type:	Mass Spectrometer	Liner Type:	4 mm ID Double Taper Liner
Carrier Gas:	He, 7.0 psi	Liner Part Number:	092018
Carrier Gas Flow:	1.0 mL/min.		
Constant Flow:	On		



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)


## ANALYSIS OF ROSEMARY OIL ON BPX5

**Column Part No.:** 054101

Phase: BPX5, 0.25 µm film  
 Column: 30 m x 0.25 mm ID  
 Initial Temp.: 40 °C, 1 min.  
 Rate 1: 5 °C/min to 260 °C,  
 Final Temp: 260 °C  
 Detector Type: Mass Spectrometer  
 Carrier Gas: He, 7.0 psi  
 Carrier Gas Flow: 1.0 mL/min.  
 Constant Flow: On

Average Linear Velocity: 36 cm/sec at 40 °C

Injection Mode: Split  
 Split Ratio: 200:1  
 Purge on (Split)  
 Vent Flow: 200 mL/min.  
 Injection Volume: 0.2 µL  
 Injection Temp.: 250 °C  
 Liner Type: 4 mm ID Double  
 Taper Liner  
 Liner Part Number: 092018



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)


## ANALYSIS OF LAVENDER OIL ON BPX5

**Column Part No.:** 054101

Phase: BPX5, 0.2 5 µm film  
 Column: 30 m x 0.25 mm ID  
 Initial Temp.: 40 °C, 1 min.  
 Rate 1: 5 °C/min to 260 °C,  
 Final Temp: 260 °C  
 Detector Type: Mass Spectrometer  
 Carrier Gas: He, 7.0 psi  
 Carrier Gas Flow: 1.0 mL/min.  
 Constant Flow: On

Average Linear Velocity: 36 cm/sec at 40 °C

Injection Mode: Split  
 Split Ratio: 200:1  
 Purge on (Split)  
 Vent Flow: 200 mL/min.  
 Injection Volume: 0.2 µL  
 Injection Temp.: 250 °C  
 Liner Type: 4 mm ID Double  
 Taper Liner  
 Liner Part Number: 092018




For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

## ANALYSIS OF TEATREE OIL ON BPX5

**Column Part No.:** 054101

Phase: BPX5, 0.25 µm film  
 Column: 30 m x 0.25 mm ID  
 Initial Temp.: 40 °C, 1 min.  
 Rate 1: 5 °C/min to 200 °C,  
 Final Temp: 200 °C  
 Detector Type: Mass Spectrometer  
 Carrier Gas: He, 7.0 psi  
 Carrier Gas Flow: 1.0 mL/min.  
 Constant Flow: On

Average Linear Velocity: 36 cm/sec at 40 °C  
 Injection Mode: Split  
 Split Ratio: 200:1  
 Purge on (Split)  
 Vent Flow: 200 mL/min.  
 Injection Volume: 0.2 µL  
 Injection Temp.: 250 °C  
 Liner Type: 4 mm ID Double  
 Taper Liner  
 Liner Part Number: 092018




For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

## ANALYSIS OF PINE OIL ON BPX5

**Column Part No.:** 054101  
Phase: BPX5, 0.25 µm film  
Column: 30 m x 0.25 mm ID  
Initial Temp.: 40 °C, 1 min.  
Rate 1: 5 °C/min to 260 °C,  
Final Temp: 260 °C  
Detector Type: Mass Spectrometer  
Carrier Gas: He, 7.0 psi  
Carrier Gas Flow: 1.0 mL/min.  
Constant Flow: On


Average Linear Velocity: 36 cm/sec at 40 °C  
Injection Mode: Split  
Split Ratio: 200:1  
Purge on (Split)  
Vent Flow: 200 mL/min.  
Injection Volume: 0.2 µL  
Injection Temp.: 250 °C  
Liner Type: 4 mm ID Double  
Taper Liner  
Liner Part Number: 092018



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

## ANALYSIS OF EUCALYPTUS OIL ON BPX5


<b>Column Part No.:</b>	<b>054101</b>	Average Linear Velocity:	36 cm/sec at 40 °C
Phase:	BPX5, 0.25 µm film	Injection Mode:	Split
Column:	30 m x 0.25 mm ID	Split Ratio:	200:1
Initial Temp.:	40 °C, 1 min.	Purge on (Split)	
Rate 1:	5 °C/min to 200 °C,	Vent Flow:	200 mL/min.
Final Temp:	260 °C	Injection Volume:	0.2 µL
Detector Type:	Mass Spectrometer	Injection Temp.:	250 °C
Carrier Gas:	He, 7.0 psi	Liner Type:	4 mm ID Double Taper Liner
Carrier Gas Flow:	1.0 mL/min.	Liner Part Number:	092018
Constant Flow:	On		



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# ANALYSIS OF YLANG YLANG OIL ON SOLGEL-WAX™

<b>Column Part No.:</b>	<b>054796</b>	Constant Flow:	On
Phase:	SolGel-Wax™, 0.25 µm	Average Linear Velocity:	35 cm/sec at 40 °C
Sample:	Ylang Ylang oil neat.	Injection Mode:	Split
Column:	30 m x 0.25 mm ID	Split Ratio:	120:1
Initial Temp.:	40 °C, 2 min.	Injection Volume:	0.1 µL
Rate 1:	3 °C/min to 250 °C,	Injection Temp.:	250 °C
Final Temp:	250 °C, 10 min.	Autosampler:	No
Detector Type:	Mass Spectrometer	Liner Type:	4 mm ID Double Taper Liner
Carrier Gas:	He, 25.7 psi	Liner Part Number:	092018
Carrier Gas Flow:	1.8 mL/min.	Full Scan / SIM:	Full scan 45-450



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)




# ANALYSIS OF EUCALYPTUS OIL ON SOLGEL-WAX™

**Column Part No.: 054796**

Phase: SolGel-Wax™, 0.25 µm film  
 Sample: neat  
 Column: 30 m x 0.25 mm ID  
 Initial Temp.: 40 °C, 1 min.  
 Rate 1: 8 °C/min to 220 °C,  
 Final Temp: 220 °C, 5 min.  
 Detector Type: Mass Spectrometer  
 Carrier Gas: He, 25.7 psi  
 Carrier Gas Flow: 1.8 mL/min.

Constant Flow: On  
 Average Linear Velocity: 35 cm/sec at 40 °C  
 Injection Mode: Split  
 Split Ratio: 100:1  
 Injection Volume: 0.2 µL  
 Injection Temp.: 250 °C  
 Liner Type: 4 mm ID Single Taper Liner  
 Liner Part Number: 092017  
 Full Scan / SIM: Full scan 45-450




For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# ANALYSIS OF NUTMEG OIL ON SOLGEL-WAX™

**Column Part No.:** 054796

Phase: SolGel-Wax™, 0.25 µm  
Sample: 2000-3000 ppm in ethanol  
Column: 30 m x 0.25 mm ID  
Initial Temp.: 40 °C, 1 min.  
Rate 1: 8 °C/min to 220 °C,  
Final Temp: 220 °C, 5 min.  
Detector Type: Mass Spectrometer  
Carrier Gas: He, 25.7 psi  
Carrier Gas Flow: 1.8 mL/min.

Constant Flow: On  
Average Linear Velocity: 35 cm/sec at 40 °C  
Injection Mode: Split  
Split Ratio: 100:1  
Injection Volume: 0.2 µL  
Injection Temp.: 250 °C  
Liner Type: 4 mm ID  
Single Taper Liner  
Liner Part Number: 092017  
Full Scan / SIM: Full scan 45-450

**Components**

1. α-Pinene
2. α-Phellandrene
3. β-Pinene
4. Sabinene
5. Myrcene
6. Limonene
7. γ-Terpinene
8. p-Cymene
9. Terpinen-4-ol
10. Safrole
11. Myristicin


For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# ANALYSIS OF SPEARMINT OIL ON SOLGEL-WAX™

**Column Part No.:** 054796

Phase: SolGel-WAX™, 0.25 µm film  
 Sample: neat  
 Column: 30 m x 0.25 mm ID  
 Initial Temp.: 40 °C, 1 min.  
 Rate 1: 8 °C/min to 220 °C,  
 Final Temp: 220 °C, 5 min.  
 Detector Type: Mass Spectrometer  
 Carrier Gas: He, 25.7 psi  
 Carrier Gas Flow: 1.8 mL/min.

Constant Flow: On  
 Average Linear Velocity: 35 cm/sec at 40 °C  
 Injection Mode: Split  
 Split Ratio: 100:1  
 Injection Volume: 0.2 µL  
 Injection Temp.: 250 °C  
 Liner Type: 4 mm ID  
 Single Taper Liner  
 Liner Part Number: 092017  
 Full Scan / SIM: Full scan 45-450



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# ANALYSIS OF TASMANIAN LAVENDER OIL SOLGEL-WAX™

**Column Part No.:** 054796

Phase: SolGel-WAX™, 0.25 µm film

Sample: neat

Column: 30 m x 0.25 mm ID

Initial Temp.: 40 °C, 1 min.

Rate 1: 8 °C/min to 220 °C,

Final Temp: 220 °C, 5 min.

Detector Type: Mass Spectrometer

Carrier Gas: He, 25.7 psi

Carrier Gas Flow: 1.8 mL/min.

Constant Flow: On

Average Linear Velocity: 35 cm/sec at 40 °C

Injection Mode: Split

Split Ratio: 100:1

Injection Volume: 0.2 µL


Injection Temp.: 250 °C

Liner Type: 4 mm ID

Single Taper Liner

Liner Part Number: 092017


Full Scan / SIM: Full scan 45-450



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# ANALYSIS OF TEATREE OIL ON SOLGEL-WAX™

<b>Column Part No.:</b>	<b>054796</b>	Constant Flow:	On
Phase:	SolGel-WAX, 0.25 µm film	Average Linear Velocity:	35 cm/sec at 40 °C
Sample:	10 ppm in methanol	Injection Mode:	Split
Column:	30 m x 0.25 mm ID	Split Ratio:	100:1
Initial Temp.:	40 °C, 1 min.	Injection Volume:	0.2 µL
Rate 1:	8 °C/min to 220 °C,	Injection Temp.:	250 °C
Final Temp:	220 °C, 5 min.	Liner Type:	4 mm ID
Detector Type:	Mass Spectrometer	Liner Part Number:	092017
Carrier Gas:	He, 25.7 psi	Full Scan / SIM:	Full scan 45-450
Carrier Gas Flow:	1.8 mL/min.		



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# Physicochemical Characterization of European Blackcurrants on SolGel-WAX™

## INTRODUCTION

The yield and relative distribution of flavor and fragrance components in foodstuffs, beverages and other agricultural produce varies with plant or animal species, agricultural conditions, geographic range, and climate. Production and processing techniques can also have a major effect. In the case of higher value products, especially those that are specific to a particular region, a premium price may be obtained based on the expectation that the product possesses certain gustatory or olfactory characteristics. Indeed in Europe today it is possible to gain legal protection similar to a trademark, for such regional "specialities". Consistency is everything to the manufacturers of such products, and traditionally they have often relied on a highly experienced "nose" to ensure it. Gas chromatography is also a valuable and widely used tool, capable of producing high quality data where required for commercial, or legal, reasons.

Complete separation of closely related flavour and fragrance compounds is not always possible using non-polar phases such as 100 % dimethylpolysiloxane and 5% phenyl 95 % dimethylpolysiloxane. Complementary separation based on a H-bond accepting phase is particularly useful for resolving alcohols and other polar components. Polyethylene glycol (or 'wax') phases provide the ideal polarity but are prone to oxidative phase loss, which gives poor retention time reproducibility over time, making comparisons between different samples difficult. SGE's SolGel-WAX™ columns are much more resistant to oxidative phase loss than conventional wax phases, having lower bleed and greater retention time reproducibility with use. The combination of low bleed and resistance to oxygen, make SGE's SolGel-WAX™ columns an excellent choice for the QC analysis of flavor and fragrance compounds.

This study, using GC-MS, demonstrates the characterisation of six samples of blackcurrants, from several regions across Europe.

## EXPERIMENTAL CONDITIONS

The six blackcurrant samples were analysed within 2 days of their receipt, and were refrigerated prior to extraction. Each sample was then sealed in a bottle using a stopper with a septum and warmed to 60 °C prior to extraction. SPME extraction was carried out with a Divinylbenzene/Carboxen/Poly-dimethylsiloxane StableFlex 50/30 µm fiber (Supelco p/n: 57328-U) for 15 minutes.

The desorption of the volatile components was carried out directly in the GC injector (250 °C) for 90 seconds.

SolGel-WAX™ 30 m x 0.25 mm ID x 1 µm (Part No. 054787).

## RESULTS

A typical chromatographic profile of one of the various batches of blackcurrants is presented in Figure 1.

The qualitative analysis shows the presence of 60 components in the six extracts of blackcurrant buds. The majority of these components have been previously quoted in the literature, except the 1-Hydroxy-2-methyl ethyl ketone and hydroxyphenylethyl alcohol, which are identified in the present study. Only one component, named X, could not be identified.

Of the 60 chemical components present in the six extracts, 16 components were selected to demonstrate the different characteristics of the extracts obtained. These specific components were chosen, either because they were indicative of the balance between woody aromas and the more desirable fruity ones, or because there were notable differences in their concentrations from one extract to another. Alpha-Terpinene, Benzaldehyde, Copaene, (E)-2-Caren-4-ol and Perillyl alcohol, were deleted from this list despite their specificity, because of their low-level concentrations. The 16 components selected are presented in Table 1.




Figure 1. Chromatographic profile of blackcurrant buds from Batch 1

## APPLICATION NOTE

Chemical Compounds	Retention time (min)	% of total area					
		Batch1	Batch2	Batch3	Batch4	Batch5	Batch6
3-Carene	5.4	1,63	2,26	2,04	1,00	2,78	2,48
Limonene	6.55	0,16	0,78	0,13	0,19	0,34	0,35
β-Phellandrene	6.73	0,49	2,58	2,73	0,74	1,39	1,28
4-Carene	8.72	4,98	3,26	7,04	2,59	3,01	3,02
X	14.00	2,41	1,19	0,19	2,39	3,77	3,81
Bornyl acetate	17.59	1,10	7,31	23,72	1,27	0,99	1,16
beta Caryophyllene	18.44	26,72	21,14	14,43	18,16	19,02	18,59
gamma Elemene	19.34	-	-	-	5,59	8,9	9,01
Allo-aromadendrene	19.48	1,78	1,16	3,25	0,8	-	-
alpha Caryophyllene	20.26	13,54	10,33	4,04	8,97	8,16	8,02
Germacrene D	21.23	5,17	0,64	0,67	3,98	2,36	2,28
delta-Cadinene	22.51	1,51	0,63	0,66	1,82	1,27	1,29
p-Cymen-8-ol	24.68	1,27	1,61	0,63	0,45	0,45	0,45
2-Methoxyphenol	24.82	0,12	0,29	0,11	1,88	1,74	1,78
Caryophyllene oxide	28.09	0,65	2,08	4,33	2,06	1,71	1,82
Spathulenol	31.51	0,48	1,27	2,74	3,62	2,47	2,53

Table 1: List of components chosen to distinguish between blackcurrant batches

### CONCLUSIONS

After extraction by SPME and analysis by GC-MS, using SGE's SolGel-WAX™ GC columns of the six batches, it was possible to identify 60 components, including 16 that differed from batch to batch sufficiently to allow determination of each sample's source and/or quality.

### REFERENCES

We thank M. Bernet and M. Ditsch of the ISIPCA Group, Research and Studies Centre for Fragrance, Cosmetics and Food Flavours, France, for developing this application, and also SOCOFRUIT, les COTEAUX BOURGUIGNONS, and Mr Krausz of the CEPPARM for providing the blackcurrants.

For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

#### AUSTRALIA & PACIFIC REGION

SGE Analytical Science Pty Ltd  
Toll Free: 1800 800 167  
Tel: +61 (0) 3 9837 4200  
Fax: +61 (0) 3 9874 5672  
Email: [support@sge.com](mailto:support@sge.com)

#### CHINA

SGE Shanghai Representative Office  
Tel: +86 21 6407 9382  
Fax: +86 21 6407 9386  
Email: [china@sge.com](mailto:china@sge.com)

#### MIDDLE EAST

SGE Gulf  
Tel: +971 6 557 3341  
Fax: +971 6 557 3541  
Email: [gulfsupport@sge.com](mailto:gulfsupport@sge.com)

#### EUROPE

SGE Europe Ltd  
European Head Office  
Toll Free: 00800 2790 8999  
Toll Free Fax: 00800 2626 2609  
Tel: +44 1908 568 844  
Fax: +44 1908 566 790  
Tel France: +33 1 69 29 80 90  
Fax France: +33 1 69 29 09 25  
Tel Germany: +49 (0) 6155 / 60746 0  
Fax Germany: +49 (0) 6155 / 60746 50  
Email: [europe@sge.com](mailto:europe@sge.com)

#### INDIA

SGE Laboratory Accessories Pvt Ltd  
Tel: +91 22 24715896  
Fax: +91 22 24716592  
Email: [sgeindia@vsnl.com](mailto:sgeindia@vsnl.com)

#### UNITED STATES OF AMERICA

SGE Incorporated  
Toll Free: (800) 945 6154  
Tel: +1 512 837 7190  
Fax: +1 512 836 9159  
Email: [usa@sge.com](mailto:usa@sge.com)

#### JAPAN

SGE Japan Inc  
Tel: +81 45 222 2885  
Fax: +81 45 222 2887  
Email: [japan@sge.com](mailto:japan@sge.com)



# Analysis of Serai Essential Oil by Fast GCMS

## INTRODUCTION

Analysis of essential oil by GCMS normally takes a long time due to its complex compositions. A basic analysis often takes between 60 minutes to 100 minutes for a complete identification of its compounds. The time can be shortened by carrying out the same analysis on a Fast GCMS. This application note demonstrates the analysis of west indian lemongrass, otherwise known as Serai oil, by Fast GCMS.

## Instrument Parameters (Fast GCMS):

Instrument: Shimadzu GCMS 2010 ULTRA  
 Column: BPX5 10 m x 0.10 mm ID x 0.10 µm  
 Column Part Number : 054099  
 Injection Mode: High Pressure injection (550 kPa, Time: 1.0 min) (Split 1:50)  
 Flow Control Mode: Linear Velocity (50 cm/sec)  
 Injection Temperature: 250 °C  
 Column Temperature: Temperature Gradient [Initial 50 °C (hold 1 min), 75 °C/min until 150 °C (hold 0.30 min), 60 °C/min until 200 °C (hold 0.5 min), 750 °C/min until 300 °C].  
 Interface Temperature: 300 °C  
 Ion Source Temperature: 200 °C  
 Mass Range: 40 - 500.00 m/z.  
 Scan Speed: 20,000 u/sec.

## RESULTS




Figure 1. Total Ion Chromatogram of Fast GC-MS and its compound names.

## REFERENCES

Analysis of Essential Oil by Fast GCMS by Fisher Scientific Malaysia, Ang May Yen, Chow Tan Wei, Hoh Lian Hoi, Hazni Hashim.

For more information contact our technical customer support team on: techsupport@sge.com

### AUSTRALIA & PACIFIC REGION

SGE Analytical Science Pty Ltd  
 Toll Free: 1800 800 167  
 Tel: +61 (0) 3 9837 4200  
 Fax: +61 (0) 3 9874 5672  
 Email: support@sge.com

### CHINA

SGE Shanghai Representative Office  
 Tel: +86 21 6407 9382  
 Fax: +86 21 6407 9386  
 Email: china@sge.com

### MIDDLE EAST

SGE Gulf  
 Tel: +971 6 557 3341  
 Fax: +971 6 557 3541  
 Email: gulfsupport@sge.com

### EUROPE

SGE Europe Ltd  
 European Head Office  
 Toll Free: 00800 2790 8999  
 Toll Free Fax: 00800 2626 2609  
 Tel: +44 1908 568 844  
 Fax: +44 1908 566 790  
 Tel France: +33 1 69 29 80 90  
 Fax France: +33 1 69 29 09 25  
 Tel Germany: +49 (0) 6155 / 60746 0  
 Fax Germany: +49 (0) 6155 / 60746 50  
 Email: europe@sge.com

### INDIA

SGE Laboratory Accessories Pvt Ltd  
 Tel: +91 22 24715896  
 Fax: +91 22 24716592  
 Email: sgeindia@vsnl.com

### UNITED STATES OF AMERICA

SGE Incorporated  
 Toll Free: (800) 945 6154  
 Tel: +1 512 837 7190  
 Fax: +1 512 836 9159  
 Email: usa@sge.com

### JAPAN

SGE Japan Inc  
 Tel: +81 45 222 2885  
 Fax: +81 45 222 2887  
 Email: japan@sge.com

 **SGE Analytical Science**  
 www.sge.com