


# ANALYSIS OF ACID/NEUTRAL DRUGS ON BPX35

<b>Column Part No.:</b>	<b>054711</b>	Final Temp.:	300 °C, 5 min
Phase:	BPX35, 0.25 µm	Carrier Gas:	He, 150 kpa
Column:	25 m x 0.22 mm ID	Injection Mode:	Split, (20:1)
Initial Temp.:	100 °C, 1 min	Detector:	FID, 380 °C
Rate:	10 °C/min		


Note: BPX35 is a low bleed column with a maximum temperature of 360 °C. Very compatible with GC/MS systems..



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# ANALYSIS OF BASIC DRUGS ON BPX35

<b>Column Part No.:</b>	<b>054711</b>	Final Temp.:	325 °C, 5 min
Phase:	BPX35, 0.25 µm	Carrier Gas:	Helium 150 kpa
Column:	25 m x 0.22 mm ID	Injection Mode:	Split, 0.5 µL (20:1)
Initial Temp.:	100 °C	Detector:	FID, 380 °C
Rate:	5 °C/min		




For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

## ANALYSIS OF ALKALOIDS ON BP5

<b>Column Part No.:</b>	<b>054198</b>	Final Temp:	300 °C, 0 min
Phase:	BP5, 1.0 µm film	Detector:	FID
Column:	25 m x 0.53 mm ID	Sensitivity:	128 x 10 <sup>-12</sup> AFS
Initial Temp.:	200 °C, 0 min	Injection Mode:	Split
Rate:	25 °C/min		


Note: A 0.53 mm ID column can be used to screen samples rapidly.



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# ANALYSIS OF SEDATIVES/HYPNOTICS ON BP1

<b>Column Part No.:</b>	<b>054087</b>	Final Temp.:	250 °C, 3 min
Phase:	BP1, 1.0 µm film	Detector:	FID
Column:	25 m x 0.53 mm ID	Sensitivity:	1024 x 10 <sup>-12</sup> AFS
Initial Temp.:	180 °C, 0 min	Injection Mode:	Split
Rate:	10 °C/min		



- Components**
1. Allobarbital
  2. Aprobarbital
  3. Diphenhydramine
  4. Mephobarbital
  5. Methapyrilene
  6. Chloropheniramine


For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# ANALYSIS OF BASIC DRUG SCREEN ON BPX5

(10-20 ng/component)

<b>Column Part No.:</b>	<b>054131</b>	Final Temp.:	310 °C
Phase:	BPX5, 1.0 µm	Detector:	FID
Column:	25 m x 0.53 mm I.D.	Injector:	Split, 240 °C
Initial Temp.:	120 °C	Carrier Gas:	H <sub>2</sub> , 2 psi
Rate:	10 °C/min		


Note: The low bleed nature of the BPX5 allows trace analysis to be performed.



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# STEROID ANALYSIS (UNDERIVATISED) ON BPX5

<b>Column Part No.:</b>	<b>054113</b>	Final Temp.:	350 °C, 10 min
Phase:	BPX5, 0.25 µm	Detector:	FID
Column:	25 m x 0.22 mm ID	Sensitivity:	32 x 10 <sup>-12</sup> AFS
Initial Temp.:	180 °C	Injection Mode:	Split
Rate:	8 °C/min	Carrier Gas:	H <sub>2</sub> , 10 psi




For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# ANALYSIS OF TRICYCLIC ANTIDEPRESSANTS ON BPX35

<b>Column Part No.:</b>	<b>054711</b>	Final Temp:	280 °C
Phase:	BPX35, 0.25 µm	Carrier Gas:	Helium, 150 kpa
Column:	25 m x 0.22 mm ID	Injection Mode:	Split (20:1)
Initial Temp.:	210 °C, 1 min	Detector:	FID, 380 °C
Rate:	5 °C/min		


Note: BPX35 is a low bleed, chemically inert phase which allows trace analysis to occur.



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# UNDERIVATISED BARBITURATES ON BP5

<b>Column Part No.:</b>	<b>054197</b>	Carrier Gas:	Hydrogen
Phase:	BP5, 1.0 $\mu$ m	Carrier Flow:	10 mL/min
Column:	12 m x 0.53 mm I.D.	Injection Volume:	0.1 $\mu$ L
Temp:	195 $^{\circ}$ C		




For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)



## ANALYSIS OF BLOOD ALCOHOL ON BP20

<b>Column Part No.:</b>	<b>054442</b>	Detector:	FID
Phase:	BP20, 1.0 $\mu\text{m}$ film	Sensitivity:	$64 \times 10^{-12}$ AFS
Column:	25 m x 0.32 mm ID	Injection Mode:	Split
Initial Temp:	Isothermal, 60 $^{\circ}\text{C}$		


*Note: The BP20 column allows the use of aqueous solutions.*



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# ANALYSIS OF DRUGS OF ABUSE ON BPX35


<b>Column Part No.:</b>	<b>054711</b>	Temp 2:	200 °C
Phase:	BPX35, 0.25 µm film	Rate 2:	7 °C/min
Column:	25 m x 0.22 mm ID	Temp 3:	295 °C
Initial Temp.:	80 °C	Rate 3:	20 °C/min
Rate 1:	15 °C/min	Final Temp.:	340 °C, 6 min



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

## ANALYSIS OF HORSE RACING TEST MIX ON BPX5


<b>Column Part No.:</b>	<b>054101</b>	Constant Flow:	On
Phase:	BPX5, 0.25 µm film	Average Linear Velocity:	45 cm/sec at 75 °C
Column:	30 m x 0.25 mm ID	Injection Mode:	Splitless
Horse Racing standard*:	10 ppm in methanol	Purge on Time:	0.5 min
Initial Temp:	75 °C, 2 min	Purge on (Split)	
Rate 1:	15 °C/min to 300 °C	Vent Flow:	60 mL/min
Rate 2:	20 °C/min to 320 °C	Injection Volume:	1 µL
Final Temp:	320 °C, 8 min.	Injection Temperature:	250 °C
Detector Type:	Mass Spectrometer	Liner Type:	4 mm ID
Carrier Gas:	He, 14.5 psi		Double Taper Liner
Carrier Gas Flow:	1.5 mL/min	Liner Part Number:	092018



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# ANALYSIS OF HORSE RACING TEST MIX ON SOLGEL-1ms™


<b>Column Part No.:</b>	<b>054795</b>	Carrier Gas Flow:	1.7 mL/min
Phase:	SolGel-1ms™	Constant Flow:	On
	0.25 µm film	Average Linear Velocity:	35 cm/sec at 75 °C
Column:	30 m x 0.25 mm ID	Injection Mode:	Splitless
Horse Racing standard*:	10 ppm in methanol	Purge on Time:	0.5 min
Initial Temp.:	75 °C, 2 min	Purge on (Split) Vent Flow:	60 mL/min
Rate 1:	15 °C/min to 300 °C	Injection Volume:	1 µL
Rate 2:	20 °C/min to 320 °C	Injection Temperature:	250 °C
Final Temp.:	320 °C, 8 min.	Liner Type:	4 mm ID
Detector Type:	Mass Spectrometer	Liner Part Number:	092017
Carrier Gas:	He, 28.7 psi		



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# ANALYSIS OF HORSE RACING TEST MIX ON BPX35


<b>Column Part No.:</b>	<b>054701</b>	Constant Flow:	On
Phase:	BPX35 0.25 µm film	Average Linear Velocity:	45 cm/sec at 75 °C
Column:	30 m x 0.25 mm ID	Injection Mode:	Splitless
Horse Racing standard*:	10 ppm in methanol	Purge on Time:	0.5 min
Initial Temp.:	75 °C, 2 min	Purge on (Split) Vent Flow:	60 mL/min
Rate 1:	15 °C/min to 300 °C	Injection Volume:	1 µL
Rate 2:	20 °C/min to 320 °C	Injection Temperature:	250 °C
Final Temp.:	320°C, 8 min	Liner Type:	4 mm ID Double Taper Liner
Detector Type:	MSD	Liner Part Number:	092018
Carrier Gas:	He, 14.5 psi		
Carrier Gas Flow:	1.5 mL/min		



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# ANALYSIS OF VARIOUS DRUGS ON BPX50


<b>Column Part No.:</b>	<b>054751</b>	Carrier Gas Flow:	1.8 mL/min
Phase:	BPX50, 0.25µm film	Constant Flow:	On
Column:	30 m x 0.25 mm ID	Average Linear Velocity:	35 cm/sec at 40 °C
Sample:	5-10 ppm in methanol	Injection Mode:	Splitless
Initial Temp:	150 °C, 0.5 min	Purge on Time:	0.5 min
Rate 1:	10 °C/min to 180 °C	Purge on (Split)Vent Flow:	60 mL/min
Rate 2:	1.5 °C/min to 220 °C	Injection Volume:	1 µL
Rate 2:	30 °C/min to 260 °C	Injection Temperature:	250 °C
Final Temp:	260 °C, 5 min	Liner Type:	4 mm ID Single Taper Liner
Detector Type:	FID	Liner Part Number:	092017
Detector Temp.:	320 °C	Full Scan / SIM:	Full scan 45-450
Carrier Gas:	He, 25.7 psi		



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# ANALYSIS OF A VARIETY OF ANTIDEPRESSANT AND ANTICONVULSANT DRUGS ON BPX50

<b>Column Part No.:</b>	<b>054751</b>	Carrier Gas Flow:	1.8 mL/min.
Phase:	BPX50, 0.25 µm film	Constant Flow:	On
Column:	30 m x 0.25 mm ID	Average Linear Velocity:	35 cm/sec at 40 °C
Sample:	5-10 ppm in methanol	Injection Mode:	Splitless
Initial Temp:	150 °C, 0.5 min	Purge on Time:	0.5 min
Rate 1:	10 °C/min to 180 °C	Purge on (Split) Vent Flow:	60 mL/min
Rate 2:	1.5 °C/min to 220 °C	Injection Volume:	1 µL
Rate 2:	30 °C/min to 260 °C	Injection Temperature:	250 °C
Final Temp:	260 °C, 5 min	Liner Type:	4 mm ID Single Taper Liner
Detector Type:	FID	Liner Part Number:	092017
Detector Temp.:	320 °C	Full Scan / SIM:	Full scan 45-450
Carrier Gas:	He, 25.7 psi		



For more information contact our technical customer support team on: [techsupport@sge.com](mailto:techsupport@sge.com)

# FALSE POSITIVE DRUG TESTS FROM CONTAMINATED FEED

An On-line and Off-line Application of Micro-SPE (MEPS)

## INTRODUCTION

Papaver somniferum (opium poppy) is a feed contaminant that can result in positive drug tests for racing horses.

A Micro-Extraction Packed Sorbent (MEPS) / GCMS method is described below for the regulatory testing of equine and human urine samples. The method allows for the separation of morphine and its metabolites from the metabolites of potential botanical markers that indicate the ingestion of poppy seeds or straw. MEPS is the miniaturization of conventional SPE from milliliter to microliter bed volumes that allows SPE to be used with very small samples. The manipulation of the small volumes is achieved with a precision gas tight syringe. With a typical void volume of 7  $\mu\text{L}$ , the MEPS elution is compatible with GC and LC inlets making it ideal for integration into an automated sampling system for on-line SPE.

## EXPERIMENTAL CONDITIONS

A 300  $\mu\text{L}$  sample of diluted equine urine from an animal receiving contaminated feed was hydrolyzed with  $\beta$ -glucuronidase or acid, filtered and extracted on a C8/SCX MEPS cartridge conditioned with methanol (30  $\mu\text{L}$ ), and potassium phosphate buffer (0.2 M, pH 6, 30  $\mu\text{L}$ ) at a flow rate of 5  $\mu\text{L}/\text{sec}$ . The exhausted fraction was ejected at the same rate and the sorbent washed with 100  $\mu\text{L}$  phosphate buffer, 50  $\mu\text{L}$  acetic acid (1% v/v) and 100  $\mu\text{L}$  methanol. The sorbent was dried with air (3 x 80  $\mu\text{L}$  at 50  $\mu\text{L}/\text{sec}$ ) and the sorbent eluted with 20  $\mu\text{L}$  dichloromethane-isopropanol-ammonia (49:49:2). The organic phase was evaporated under nitrogen and derivatized with 10  $\mu\text{L}$  of acetic anhydride-pyridine (1:2) at 80  $^{\circ}\text{C}$  for 30 minutes before evaporation and reconstitution in 5  $\mu\text{L}$  of ethyl acetate.

## RESULTS

The extract was analyzed by GCMS on a BPX5 column (Figure 1a and 1b).

## CONCLUSIONS

This application of mixed mode C8/SCX MEPS for a complex biological fluid, allowed the microscale preparation of a small volume sample with comparable performance to conventional SPE techniques. Used off-line with derivatization and GCMS here, the sample was also suitable for on-line ESI-LCMSMS analysis by changing the elution solvent to methanol-ammonia (98:2) or methanol-trimethylamine (98:2).

In most cases, MEPS allows the same level of sample concentration as is possible with off-line conventional SPE while providing opportunities for truly hybrid multi-dimensional methods. MEPS methods may be readily adapted from established SPE methods including those based on mixed mode or complex chemistries.


Like SPE, MEPS is for use with liquid samples (either normal or reversed phase) and yields four fractions: the unretained, weakly bound, strongly bound and irreversibly bound. However, because MEPS is a double pass system (sample and solvent enter and exit from the bottom of the bed), the weakly bound fraction (commonly the interferences eliminated by washing) is less strongly bound. The irreversibly bound fraction affects MEPS and conventional SPE and is usually associated with sorbent wetting rather than sample purification and so the irreversible binding of matrix material from one sample does not preclude reuse of the device for a sample of the same type.

Also like conventional SPE, the number of times the device can be re-used is dependent on the sample matrix. For simple applications, MEPS devices have been used successfully for more than 50 cycles.


## REFERENCES

(1) Wynne PM, The accidental contamination of animal feed by naturally occurring opiates with particular reference to morphine. An independent review and analysis of current knowledge. The British Equestrian Trade Association (London), 2005; pp1-143.





**Figure 1a:** The analysis of *Papaver somniferum* ssp. *setigerum* by GCMS following micro-SPE extraction of poppy residues and GCMS analysis of the recovered fraction.



**Figure 1b:** The analysis of enzyme hydrolyzed horse urine by GCMS following micro-SPE extraction and micro-derivatization (peracetylation) of the recovered basic fraction. (Inset m/z 367 corresponds to O-desmethylpapaverine metabolites and m/z 395 is the O,O'-didesmethylmetabolites)