

Thermally Stable Polar Stationary Phase

Introduction


- Excellent separation of fatty acid methyl ester (FAME) geometric and positional isomers
- Suitable for isomeric drug compounds and carbohydrates
- Equivalent polarity to 70% cyanopropyl siloxane
- Modified silphenylene siloxane backbone for high thermal stability
- Low bleed

BPX70 is a phase specifically designed to aid in the analysis of FAMES no matter how simple or complex the mixture. BPX70 utilizes the phase technology developed by SGE to enhance the performance of what would be, "just another cyanopropyl stationary phase".

The incorporation of the aromatic silphenylene group into the siloxane backbone provides many benefits. As the initials BP (Bonded Phase) imply, the stationary phase is

crosslinked. The most obvious benefit of crosslinking is in allowing the column to be rinsed to remove contaminant material. Crosslinking also has a stabilizing effect on the film coating which, for high cyanopropyl content phases, can undergo break down at elevated temperatures.

Bis-cyanopropylsiloxane - silphenylene polymer




Analysis of a Complex FAME mixture on BPX70

Phase: BPX70, 0.25µm film
Sample: 1000 ppm in dichloromethane
Column: 30m x 0.25mm ID
 Initial Temp.: 50°C, 2min
 Rate 1: 4°C/min to 250°C
 Final Temp.: 250°C, 15min
 Detector Type: Mass Spectrometer
 Carrier Gas: He, 7.63 psi
 Carrier Gas Flow: 0.5mL/min
 Constant Flow: On

Average Linear Velocity: 19 cm/sec at 50°C
 Injection Mode: Split
 Split Ratio: 100:1
 Injection Volume: 1µL
 Injection Temp.: 250°C
 Liner Type: 4mm ID Single Taper Liner
Liner Part No.: 092017
Column Part No.: 054622
 ms-NoVent™ Part No.: 113400
 HP5973 restrictor: 113409
 Full Scan / SIM: Full scan 45-500

1. C4:0
2. C6:0
3. C8:0
4. C10:0
5. C11:0
6. C12:0
7. C13:0
8. C14:0
9. C14:1
10. C15:0
11. C15:1
12. C16:0
13. C16:1
14. C17:0
15. C17:1
16. C18:0
17. C18:1n9 (*trans*)
18. C18:1n9 (*cis*)
19. C18:2n6 (*trans*)
20. C18:2n6 (*cis*)
21. C18:3n6
22. C18:3n3
23. C20:0
24. C20:1n9
25. C21:0
26. C20:3n6
27. C20:3n3
28. C20:4n6
29. C22:0
30. C22:1n9
31. C20:5n3
32. C23:0
33. C22:2
34. C24:0
35. C24:1n9
36. C22:6n3




However, the most dramatic improvement is the increase in the maximum operating temperature of the phase. Suddenly those analyses which were difficult due to temperature limitations, become routine. With a maximum operating temperature of 260°C, continual operation at 250°C (generally regarded as the upper limit for other polar columns) has a minimal effect to the life of a BPX70 column.

Specifications

BPX70 is a modified polycyanopropylsiloxane phase with equivalent polarity to 70% cyanopropyl.

Minimum Operating Temp.= 50°C
 Maximum Continuous Temp.= 250°C
 Maximum Cycling Temp.= 260°C

FAME reference standard on BPX70



Phases of Similar Polarity to BPX70



SGE	BPX70	Bonded (70% cyanopropyl (modified))
J&W	DB-23	Bonded (50% cyanopropyl)
Restek	Rtx-2330	Non-bonded
Chrompack	CP-Sil 88	Non-bonded
Supelco	SP-2330	Non-bonded

Application Areas

Although BPX70 phase was specifically designed for the separation of geometric and positional isomer of fatty acid methyl esters, wider applications have been found in areas such as the analysis of isomeric drug compounds, including tricyclic antidepressants and carbohydrates.

BPX70 columns are essential for laboratories involved in the following areas:


- Food/Nutrition analyzing Vegetable Oils (Peanut, Olive and Palm Oils)
- Fish Oils
- Dairy Products
- Medical Diagnostics
- Cosmetic Oils
- Soaps and Detergents
- Pharmaceutical Products
- Paints and Varnishes
- Plastics
- Geochemical Samples
- Sugar Products



RAPASEED OIL FAME

Phase: **BPX70, 0.25µm film**
 Column: **12m x 0.22mm ID**
 Temp.: Isothermal at 175°C
 Carrier Gas: He, 15psi
 Detector: FID
 Sensitivity: 32 x10-12AFS
 Injection Mode: Split, 50:1
 Part No.: **054601**


- 1. C14:0
- 2. C16:0
- 3. C18:0
- 4. C18:1
- 5. C18:2
- 6. C18:3
- 7. C20:0
- 8. C20:1
- 9. C22:0
- 10. C22:1
- 11. C24:0



FAME MIXTURE

Phase: **BPX70, 0.25µm film**
 Column: **25m x 0.32mm ID**
 Temp.: Isothermal at 170°C
 Carrier Gas: H2, 6psi
 Detector: FID
 Sensitivity: 32 x10-12AFS
 Injection Mode: Split, 50:1
 Part No.: **054606**


- 1. 16:0
- 2. 16:1 (trans)
- 3. 16:1 (cis)
- 4. 18:0
- 5. 18:1 (trans)
- 6. 18:1 (cis)
- 7. 18:2 (trans)
- 8. 18:2 (cis)
- 9. 18:3 (cis)
- 10. 20:0
- 11. 20:1 (cis)
- 12. 22:0
- 13. 22:1 (cis)
- 14. 24:0
- 15. 24:1 (cis)




BACTERIAL ACID METHYL ESTERS

Phase: **BPX70, 0.25µm film**
 Column: **25m x 0.22mm ID**
 Initial Temp.: 100°C, 1min
 Program Rate: 5°C/min
 Final Temp.: 220°C, 0min
 Carrier Gas: H2, 10psi
 Detector: FID
 Sensitivity: 32 x10-12AFS
 Injection Mode: Split, 50:1
 Part No.: **054602**

- 1. C11:0
- 2. C12:0
- 3. C13:0
- 4. C10:0 2-OH
- 5. C14:0
- 6. C15:0-i
- 7. C15:0-a
- 8. C15:0
- 9. C12:0 2-OH
- 10. C16:0-i
- 11. C16:0
- 12. C15:1n7(cis)
- 13. C12:0 3-OH
- 14. C17:0-a
- 15. C17:0
- 16. C9,10-methylene 16:0
- 17. C14:0-2-OH
- 18. C18:0
- 19. C18:1n9 (trans)
- 20. C18:1n9 (cis)
- 21. C18:1n7 (cis)
- 22. C18:2n6 (cis)
- 23. C19:0
- 24. C16:0 2-OH
- 25. C9,10-methylene 18:0
- 26. C20:0
- 27. C18:0 2-OH




Analysis of Sugar Alditol Acetates on BPX70



Phase: BPX70, 0.25µm film
Column: 12m x 0.32mm ID
Temp.: Isothermal at 210°C
Carrier Gas: H, 3.3psi 2
Detector: FID
Sensitivity: 16 x10-12AFS
Injection Mode: Split, 50:1
Part No.: 054605

1. Erythritol
2. Deoxyribitol
3. Rhamnitol
4. Fucitol
5. Ribitol
6. Arabinitol
7. Xylitol
8. Deoxyglucitol
9. Allitol
10. Mannitol
11. Galactitol
12. Glucitol
13. Inositol

FATTY ACID METHYL ESTER MIXTURE




Phase: BPX70, 0.25µm film
Column: 25m x 0.22mm ID
Initial Temp.: 100°C, 0min
Program Rate: 4°C/min
Final Temp.: 260°C, 5min
Detector: FID
Sensitivity: 128 x10-12AFS
Injection Mode: Split, 50:1
Part No.: 054602

1. 13:0
2. 16:0
3. 16:1 (*trans*)
4. 16:1 (*cis*)
5. 17:0
6. 18:0
7. 18:1 (*trans*)
8. 18:1 (*cis*)
9. 18:2 (*trans*)
10. 18:2 (*cis*)
11. 20:0
12. 20:1
13. 20:2
14. 20:3
15. 20:4
16. 25:0
17. 26:0
18. 28:0
19. 30:0
20. 31:0

Analysis of a 24 Component FAME Mixture on BPX70

Phase: BPX70, 0.25µm film
Sample: 100ppm in Hexane
Column: 30m x 0.25mm ID
 Initial Temp.: 150°C, 0.5 min
 Rate 1: 10°C/min to 180°C
 Rate 2: 1.5°C/min to 220°C
 Rate 3: 30°C/min to 260°C
 Final Temp.: 260°C, 5min
 Detector Type: FID
 Detector Temp.: 280°C
 Carrier Gas: He, 20.7 psi
 Carrier Gas Flow: 1.3 mL/min
 Constant Flow: On
 Average Linear Vel.: 35 cm/sec at 150°C
 Injection Mode: Split
 Split Ratio: 15:1
 Injection Volume: 1µL
 Injection Temp.: 250°C
 Autosampler: No
 Liner Type: 4mm ID Single Taper Liner
Liner Part No.: 092017
Column Part No.: 054622

1. C12:0
2. C14:0
3. C14:1ω5
4. C16:0
5. C16:1ω7
6. C18:0
7. C18:1ω9
8. C18:2ω6
9. C18:3ω6
10. C18:3ω3
11. C20:0
12. C20:1ω9
13. C20:2ω6
14. C20:3ω9
15. C20:3ω6
16. C20:4ω6
17. C22:0
18. C22:1ω9
19. C20:5ω3
20. Unknown
21. C22:4ω6
22. C24:0
23. C24:1ω9
24. C22:5ω3
25. C22:6ω3



ORDERING INFORMATION - BPX70

Column Type	ID (mm)	Length (m)	Film Thickness (µm)	Temperature Limits (°C)	Part No.
BPX70	0.10	10	0.20	50 to 250/260	054600
		12	0.25	50 to 250/260	054601
		25	0.25	50 to 250/260	054602
		30	0.25	50 to 250/260	054612
		50	0.25	50 to 250/260	054603
		60	0.25	50 to 250/260	054613
	0.25	15	0.25	50 to 250/260	054621
		30	0.25	50 to 250/260	054622
		60	0.25	50 to 250/260	054623
		120	0.25	50 to 250/260	054624
	0.32	12	0.25	50 to 250/260	054605
		25	0.25	50 to 250/260	054606
		30	0.25	50 to 250/260	054616
		50	0.25	50 to 250/260	054607
	0.53	60	0.25	50 to 250/260	054617
		15	0.50	50 to 250/260	054619
		25	0.50	50 to 250/260	054610
			30	0.50	50 to 250/260



www.sge.com

BPX70 was developed in conjunction with C.S.I.R.O. Australia.

SGE International Pty Ltd
 Toll Free: 1800 800 167
 Tel: +61 (0) 3 9837 4200
 Fax: +61 (0) 3 9874 5672
 email: support@sge.com

SGE Europe Ltd (UK)
 Tel: +44 (0) 1908 568844
 Fax: +44 (0) 1908 566790
 email: uk@sge.com

SGE (Italia) srl
 Tel: +39 (06) 4429 0206
 Fax: +39 (06) 4429 0724
 email: sge.italia@tin.it

SGE GmbH (Germany)
 Tel: +49 (0) 6151 860486
 Fax: +49 (0) 6151 860489
 email: germany@sge.com

SGE Laboratory Accessories Pvt Ltd (India)
 Tel: +91 (022) 471 5896
 Fax: +91 (022) 471 6592
 email: sgeindia@vsnl.com

SGE Incorporated (USA)
 Toll Free: 800 945 6154
 Tel: +1 512 837 7190
 Fax: +1 512 836 9159
 email: usa@sge.com

SGE Europe Ltd (France)
 Tel: +33 (0) 1 69 29 80 90
 Fax: +33 (0) 1 69 29 09 25
 email: france@sge.com

SGE Japan Inc
 Tel: +81 (045) 222 2885
 Fax: +81 (045) 222 2887
 email: japan@sge.com

China Service Centre
 Tel: +86 (10) 6588 8666
 Fax: +86 (10) 6588 6577

©Copyright SGE International Pty Ltd.
 All rights reserved.
 PD-0032-C Rev:05 06/04