



Precise Numerical Results

Ochratoxin A is a mycotoxin produced by the fungus *Aspergillus ochraceous* and also by several species of *Penicillium fungi*. Frequently found in green coffee, ochratoxin A will remain present even after the roasting process. Ochratoxin A (OTA) is also commonly found in wine, baby food, raw and processed cereals, and dried vine fruits—all of which are regulated by European limits as low as 0.5 ppb. A potent kidney toxin and known carcinogen, OTA finds its way into human and animal diets not only from EU-regulated commodities, but also from barley, oats, and sorghum. OchraTest WB increases antibody volume to three times the level of standard VICAM columns, facilitating faster flow through the column and more efficient HPLC or LC/MS/MS sample preparation in your laboratory.

Using monoclonal antibody-based affinity chromatography, OchraTest™ WB is a new HPLC, UPLC and LC/MS/MS test for the detection of ochratoxin A using wide bore immunoaffinity columns. With a total volume of 3 mL, compared to the total of 1 mL in VICAM's standard columns, OchraTest WB allows for a faster flow rate preferred by many laboratories and is the ideal cleanup step for any HPLC.

OchraTest WB increases antibody volume to three times the level of standard VICAM columns, facilitating faster flow through the column and more efficient HPLC or LC/MS/MS sample preparation in your laboratory.



BENEFITS

Exclusive – Specifically for HPLC, UPLC or LC/MS/MS use

Durable – Long shelf life; requires no refrigeration

Versatile – For use with a variety of samples

Quick – 10 minutes to isolate toxin*

Wide Range – Detects levels as low as .25 ppb and as high as 100 ppb

Fast Flow – Passes more volume over the column

HERE'S HOW SIMPLE OCHRATEST WB IS TO USE:

EXTRACT SAMPLE

- Grind and weigh sample
- Add salt and methanol/water mixture
- Filter

DILUTE AND FILTER

- Dilute extract
- Filter


ABSORB AND ELUTE

- Pass filtrate over OchraTest WB affinity column
- Wash column
- Elute ochratoxin A and collect in a cuvette

MEASURE

- Dilute and inject eluate into HPLC or UPLC
- Determine ochratoxin A concentration

* Excluding preparation and extraction.



ORDERING INFORMATION

Cat. No.	Description	Cat. No.	Description
G1033	OchraTest WB Columns (25 per Box)	21020	Single Position Pump Stand*
G1034	OchraTestWB Columns (50 per Box)	23040	Vortex Mixer
31240	VICAM Fluted Filter Paper, 24 cm (100)	20501	500 mL Bottle Dispenser for Methanol (0-3 mL range)
31955	Microfiber Filters, 1.5 µm, 11 cm (100)	36020	Filter Funnel, 65 mm (10 per Pack)
34000	Disposable Cuvettes (250 per Pack)	36022	Filter Funnel, 105 mm (4 per pack)
35016	Methanol, HPLC Grade (4 X 4 L)	G4033	Micro-pipettor, 1 mL
36010	Disposable Plastic Beakers (25 per pack)	20656	Micro-pipette Tips for 1 mL Micro-pipettor (100)
G1113	10X concentrate of PBS (150 mL)		
G1124	Noniodized Sodium Chloride (Salt, NaCl)		
G1118	WB Column Coupling (6)		
20050	Graduated Cylinder, 50 mL		
20250	Graduated Cylinder, 250 mL		
20100	Digital Scale with AC Adapter		
20200	Commercial Blender with Stainless Steel Container		
20700	Wash Bottle, 500 mL		
21010	Cuvette Rack		

Also required; must be sourced locally:

Distilled, Reverse Osmosis, or Deionized Water
Centrifuge

* Multiple position pump stands available

Subject to change without notice.

© 2010 Waters Corporation. Waters, The Science of What's Possible, VICAM, UPLC, and OchraTest are trademarks of Waters Corporation.

The analytical methods presented in this data sheet have been researched and developed by Waters to be used exclusively with OchraTest products. These analytical methods have been validated in the Waters laboratories to perform to the specifications indicated in the OchraTest procedures. The user assumes all risk in using OchraTest analytical procedures and products. Waters makes no warranty of any kind, expressed or implied, other than that OchraTest products conform to Waters' printed specification and quality control standards. Waters will, at its option, repair or replace any product, or part thereof, which proves to be defective in workmanship or material. VICAM's undertaking to repair or service such products is exclusive and is in lieu of all other warranties whether written, oral, expressed, or implied, including any implied warranty of merchantability or fitness for a particular purpose. Waters shall have no liability for anticipated or lost profits or any loss, inconvenience or damage whether direct, indirect, incidental, consequential or otherwise, to person or property, or for strict liability or negligence arising from or in connection with the use of these assay procedures or OchraTest products.

Orders:

1848 N Deffer Drive
Nixa, MO 65714 USA

Tel: +1.877.228.4244
+1.417.725.6588

Fax: +1.417.725.6102

orders@vicam.com

Headquarters:

34 Maple Street
Milford, MA 01757 USA

Tel: +1.800.338.4381
+1.508.482.4935

Fax: +1.508.482.4972

www.vicam.com