



Medicine beakers, mercury collector



Medicine measures, polypropylene



- Translucent and robust graduated medicine measures
- Nest to save space
- Ideal for use when taking pH measurements

Product Code	Capacity (ml)	Grads. ml	O.D. (mm)	Height mm	Pack qty.
MWY514	30	1.0	37	42	1
MWY518	50	5.0	50	53	1



Mercury collector, polyethylene



- Spilled mercury is collected by pressing the foam pad in the screw top closure against the mercury
- The mercury is drawn into and encapsulated in the pad by releasing the pressure
- The mercury is released by screwing the closure on to the jar

Product Code	Description	Pack qty.
MWY850	Mercury collector	1
MWY851	Replacement pad	1