

# Baskets

## Caleva Compatible Baskets

Part No.	Description	OEM Ref
PSBSK010-CA	10 Mesh Stainless Steel Basket	0810123
PSBSK020-CA	20 Mesh Stainless Steel Basket	0810021
PSBSK040-CA	40 Mesh Stainless Steel Basket	0810019
PSBSKSUP-CA	Suppository Basket, Plastic	–



PSBSK010-CA



PSBSK020-CA



PSBSK040-CA



PSBSKSUP-CA

## Distek Compatible Baskets

Part No.	Description	OEM Ref
PSBSK010-DK	10 Mesh Stainless Steel Basket	2821-0160
PSBSK010-DKC	10 Mesh Stainless Steel Basket, Evolution Series	2800-0041
PSBSK020-DK	20 Mesh Stainless Steel Basket	2821-0159
PSBSK020-DKC	20 Mesh Stainless Steel Basket, Evolution Series	2800-0042
PSBSK040-DK	40 Mesh Stainless Steel Basket	2821-0072
PSBSK040-DKC	40 Mesh Stainless Steel Basket, Evolution Series	2800-0032
PSBSKPTD-40	40 Mesh PTFE Coated Basket	–
PSBSKSUP-DK	Suppository Basket, Plastic	3250-0445



PSBSK010-DK



PSBSK010-DKC



PSBSK020-DK



PSBSK020-DKC



PSBSK040-DK



PSBSK040-DKC



PSBSKSUP-DK



## Erweka *Compatible Baskets*

Part No.	Description	OEM Ref
PSBSK020-EW	20 Mesh Stainless Steel Basket	90-000-0015
PSBSK040-EW	40 Mesh Stainless Steel Basket	90-000-0040
PSBSKSUP-EW	Suppository Basket, Plastic	-



PSBSK020-EW



PSBSK040-EW



PSBSKSUP-EW

## Hanson *Compatible Baskets*

Part No.	Description	OEM Ref
PSBSK010-HR	10 Mesh Stainless Steel Basket	65-220-010
PSBSK020-HR	20 Mesh Stainless Steel Basket	65-220-020
PSBSK040-HR	40 Mesh Stainless Steel Basket	65-220-000
PSBSKPTH-40	40 Mesh PTFE Coated Basket	-
PSBSKSUP-HR	Suppository Basket, Plastic	65-700-048



PSBSK010-HR



PSBSK020-HR



PSBSK040-HR



PSBSKPTH-40



PSBSKSUP-HR

All baskets are laser marked with a unique serial number for easy identification



# Baskets

## Pharmatest Compatible Baskets

Part No.	Description	OEM Ref
PSBSK010-PT	10 Mesh Stainless Steel Basket	311-2804
PSBSK020-PT	20 Mesh Stainless Steel Basket	–
PSBSK040-PT	40 Mesh Stainless Steel Basket	311-2801-6
PSBSKSUP-PT	Suppository Basket, Plastic	–



PSBSK010-PT



PSBSK020-PT



PSBSK010-SX



PSBSK020-SX



PSBSK040-PT



PSBSKSUP-PT



PSBSK040-SX



PSBSKSUP-SX

## Sotax Compatible Baskets

Part No.	Description	OEM Ref
PSBSK010-SX	10 Mesh Stainless Steel Basket	–
PSBSK020-SX	20 Mesh Stainless Steel Basket	9353
PSBSK040-SX	40 Mesh Stainless Steel Basket	2830
PSBSKSUP-SX	Suppository Basket, Plastic	–

## Zymark Compatible Baskets

Part No.	Description	OEM Ref
PSBSK040-ZM	40 Mesh Stainless Steel Basket	63127
PSBSK040-ZM2	40 Mesh Precision Length Stainless Steel Basket	–



PSBSK040-ZM



PSBSK040-ZM2

### Precision Length Basket

The Precision Length Basket (PSBSK040-ZM2) is specifically designed for use with Caliper Multidose and other robotic dissolution systems. It has a long usable life in the demanding robotic environment.



## VanKel *Compatible Baskets*

Part No.	Description	OEM Ref
PSBSK010-01	10 Mesh Stainless Steel Basket	12-2125
PSBSK020-01	20 Mesh Stainless Steel Basket	12-2120
PSBSK040-01	40 Mesh Stainless Steel Basket	12-2100
PSBSK040-01PL	40 Mesh Precision Length Stainless Steel Basket	-
PSBSKPTV-40	40 Mesh PTFE Coated Basket	12-2110
PSBSK100-01	100 Mesh Stainless Steel Basket	12-2154
PSBSK150-01	50x150 Mesh Stainless Steel Basket	12-2152
PSBSKSUP-01	Suppository Basket, Plastic	12-2130
PSSTABSK-10	Stationary Basket Assembly with Basket Shaft and 10 Mesh Basket	12-2061
PSSTABSK-20	Stationary Basket Assembly with Basket Shaft and 20 Mesh Basket	12-2060
PSSTABSK-40	Stationary Basket Assembly with Basket Shaft and 40 Mesh Basket	12-2063



PSBSK010-01



PSBSK020-01



PSBSK040-01



PSBSKPTV-40



PSBSKSUP-01



PSSTABSK-10




PSSTABSK-20



PSSTABSK-40

USP 26 Felodipine Test Basket – See page 40





## Dissolution Basket Storage

### Prolong the life of baskets and save money

To ensure they meet USP requirements, it is essential that for all tests, baskets are not bent or deformed.

The chances of accidentally damaging baskets can be significantly reduced with the use of a good quality basket holder, such as this.

#### Key Benefits:

- **Compatible** with all popular basket types (Caleva, Distek, Erweka, Hanson, Pharmatest, Sotax, VanKel, Zymark)
- O-rings ensure baskets are **held firmly** to the holder
- **Stores up to 12** baskets
- **Fully stackable**, letting you transport many baskets at a time
- **Low cost** – it pays for itself by preventing damage to just one basket!



Part No.	Description	OEM Ref
PSBSKHLD-12	Stackable Dissolution Basket Holder, 12 Position	-

# Felodipine Test Basket



## Felodipine Stationary Basket

Part No.	Description	OEM Ref
PSSTABSK-VK8	USP 26 Felodipine Test Stationary Basket	12-2069 <b>V</b>



## Compatible Vessel Covers (See page 22)

Part No.	Description	OEM Ref
PSCOVDHG-SB	Clear Stationary Basket Apparatus Hinged Cover for Distek Baths	3250-0081 <b>D</b>
PSCOVHGV-SB	Clear Stationary Basket Apparatus Hinged Cover for VanKel Baths	-

## Stationary Basket Spacer (See page 41)

Part No.	Description	OEM Ref
PSDEPSET-SB	Stationary Basket Spacer, 1cm	-



PSDEPSET-SB

**D** = Distek    **V** = VanKel