

Transdermal Testing

Transdermal Patch Holders and Accessories

Part No.	Description	OEM Ref	
PSAPPFIVE-HR	Watch Glass Transdermal Patch Holder with 17 Mesh Screen and 4 Plastic Clips	65-190-093	H
PSAPPFIVE-HRCLIP	Plastic Clip for PSAPPFIVE-HR (each)	65-190-205	H
PSAPPFIVE-HRSCR	17 Mesh PTFE Screen for PSAPPFIVE-HR	65-190-190	H
PSAPPFIVE-V35	Transdermal Patch Holder, 1.38" (35mm) Opening (40 Mesh Screen)	12-4200	V
PSAPPFIVE-V56	Transdermal Patch Holder, 2.20" (56mm) Opening (120 Mesh Screen)	12-4230	V
PSAPPFIVE-V80	Transdermal Patch Holder, 3.15" (80mm) Opening (16 x 16 Mesh Screen)	-	
PSAPP5S16-35	16 Mesh 316 Stainless Steel Screen for PSAPPFIVE-V35 (pack of 25)	-	
PSAPP5S20-35	20 Mesh 316 Stainless Steel Screen for PSAPPFIVE-V35 (pack of 25)	-	
PSAPP5S30-35	30 Mesh 316 Stainless Steel Screen for PSAPPFIVE-V35 (pack of 25)	-	
PSAPP5S40-35	40 Mesh 316 Stainless Steel Screen for PSAPPFIVE-V35 (pack of 25)	-	



PSAPPFIVE-HR



PSAPPFIVE-HRCLIP



PSAPPFIVE-HRSCR



PSAPPFIVE-V80



PSAPPFIVE-V56



PSAPPFIVE-V35



PSAPP5S16-35



PSAPP5S20-35



PSAPP5S30-35



PSAPP5S40-35



H = Hanson **V** = VanKel



PSAPPSIX-01



PSAPPSIX-EX



Rotating Cylinder
viewed from below

Rotating Cylinder

Part No.	Description	OEM Ref	
PSAPPSIX-01	Rotating Cylinder for Transdermal Patches, Electropolished 316 Stainless Steel	12-1360 65-590-011	V H
PSAPPSIX-EX	Rotating Cylinder Extension for PSAPPSIX-01, Electropolished 316 Stainless Steel	- 65-590-012	H